

Pelco BioWave®Pro+

Microwave Processor containing **Pelco ColdSpot®Pro** cooling system

Applications

- Sample preparations for EM (fixation, contrastation, dehydration, embedding, polymerisation)

Samples

- histological tissues
- suspension cells
- hard-to infiltrate biological specimens

Specifications

- Adjustable timing
- Adjustable wattages (0 – 750 W)
- Adjustable temperature (0 – 100°C)
- Vacuum On/Off system
- Bubbler On/Off system
- Stirrer On/Off system
- Temperature check probe
- Manual or Automatic running mode
- Data recording
- Protocols storage